ACCURACY AND PRECISION STUDY OF
UNISTRIP1™ TEST STRIPS
AND
OneTouch® Ultra® Test Strips
with the
OneTouch® Ultra®, Ultra®2, UltraMini® and
UltraSmart®
Blood Glucose Meters

Study conducted by Neng-Chun Diabetes Clinical, Taiwan, ROC

INTRODUCTION

Diabetes has reached epidemic conditions in the USA. According to the American Diabetes Association, 25.8 million Americans had diabetes in 2011, which is 8.3% of the population. In 2007 diabetes was listed either as the cause or contributing factor for over 230,000 deaths. Blood glucose monitoring by people with diabetes has been proven to be effective in controlling glucose values. Among the benefits of monitoring is that patients can understand the symptoms of hyper and hypoglycemia, and the impact of diet, activity and medication on glucose levels. Most people with diabetes can achieve a high level of control on their blood glucose levels by following a physician's plan that includes regular monitoring.

Most manufacturers of blood glucose monitors desire to increase the level of compliance through testing. Improvements have been made, including the introduction of no-coding systems, audible meters, faster processing times, and smaller blood sample requirements.

Though systems are better today, and more people with diabetes are aware of the benefits of testing, the high cost of name brand monitor systems continues to have a significant impact and negative deterrent to testing. The US FDA has recently cleared the UniStrip1™ Test Strips for use with the LifeScan OneTouch® Ultra®,

Ultra®2, UltraMini® and UltraSmart® meters. With "generic" test strips now available for use by people who own the OneTouch® Ultra® meters, people with diabetes have a choice to purchase the much lower priced UniStrip1™ Test Strips.

This study was designed to gather evidence regarding the performance of the OneTouch® Ultra® meters with the OneTouch® Ultra® test strips, and the performance of lower cost UniStrip1™ Test Strips that are intended and cleared by the US FDA for use with the meters.


SITE AND SUBJECTS

This clinical study was conducted at one site.

Three lots of UniStrip1™ Test Strips and one lot of OneTouch® Ultra® Test Strips were used in this evaluation. Quality control was performed on all meter systems per manufacturer-recommended instructions. The Yellow Springs Instrument (YSI) glucose reference analyzer's calibration was checked on a daily basis, with standards provided by the manufacturer. Duplicate samples were obtained for the YSI and each test strip tested.

The data collection sequence was as follows: UniStrip1™ with OneTouch® Ultra®, Ultra®2, UltraMini® and UltraSmart® Blood Glucose Meters, YSI and HCT for 100 subjects. OneTouch® Ultra® Blood Glucose Monitoring System, YSI and HCT for 52 subjects. Blood was obtained from subjects' fingertips only.

Data Collection Sequence


The subjects were mostly people with type 2 diabetes, as well as people with type 1 diabetes, and ranged in age from 25 to 78 years old. There were slightly more females than males in this study. A broad range of glucose values were obtained for this evaluation. Subjects' hematocrit ranges were from 36%-52%.


Gender		
ONETOUGH® Ultra®	Male	44.2 %
	Female	55.8 %
UniStrip1™	Male	44.0 %
	Female	56.0 %

Diabetes Type		
ONETOUGH® Ultra®	Type 1	25.0 %
	Type 2	75.0 %
UniStrip1™	Type 1	29.3 %
	Type 2	70.7 %

Age distribution of the subjects enrolled in study






Glucose distribution of the subjects enrolled in the study

RESULTS USING PARKES ERROR GRID

The UniStrip1™ with OneTouch® Ultra®, Ultra®2, UltraMini® and UltraSmart® Blood Glucose Meters provided accurate results, that are closely comparable to the results obtained using the four OneTouch® Ultra® Blood Glucose Monitoring Systems. Each system's performance data has been plotted on a Consensus Error Grid for comparison.

Strip / Meter System	% Result in zone A	% Result in zone B
UniStrip1™ (Ultra®)	97.0	3.0
UniStrip1™ (Ultra®2)	98.5	1.5
UniStrip1™ (UltraMini®)	98.0	2.0
UniStrip1™ (UltraSmart®)	100	0.0
ONETOUCH® Ultra®	96.0	4.0

RESULTS USING ISO GUIDELINES

To further measure meter performance, all systems were evaluated by using the International Standards Organization (ISO) 15197:2013 criteria for measuring accuracy in blood glucose monitoring devices. ISO combines two sets of accuracy criteria based on glucose value; values less than 100 mg/dL and values greater than or equal to 100 mg/dL.

Summary of all systems' accuracy against the YSI reference method **for YSI readings below 100 mg/dL (calculated plasma readings).**

Strip / Meter System	Within ±15 mg/dL
UniStrip1™ (Ultra®)	34/34 (100 %)
UniStrip1™ (Ultra®2)	34/34 (100 %)
UniStrip1™ (UltraMini®)	33/34 (97 %)
UniStrip1™ (UltraSmart®)	34/34 (100 %)
ONETOUCH® Ultra®	27/27 (100 %)

Summary of all systems' accuracy against the YSI reference method **for YSI readings above 100 mg/dL (calculated plasma readings).**

Strip / Meter System	Within ±15 %
UniStrip1™ (Ultra®)	64/66 (97.0 %)
UniStrip1™ (Ultra®2)	65/66 (98.5 %)
UniStrip1™ (UltraMini®)	65/66 (98.5 %)
UniStrip1™ (UltraSmart®)	66/66 (100 %)
ONETOUCH® Ultra®	24/25 (96.0 %)

Precision was evaluated by the standard deviation calculated using the difference between duplicate samples. The UniStrip1™ with OneTouch® Ultra®, Ultra®2, UltraMini® and UltraSmart® Blood Glucose Meters has a similar average deviation to the OneTouch® Ultra® Blood Glucose Monitoring System when using an Ultra® test strip.

Strip / Meter System	Average deviation (mg/dL)
UniStrip1™ (Ultra®)	1.5
UniStrip1™ (Ultra®2)	3.4
UniStrip1™ (UltraMini®)	1.5
UniStrip1™ (UltraSmart®)	3.3
ONETOUCH® Ultra®	3.5

Note: The standard deviation of duplicates is an assessment of the system's reproducibility. The formula used to calculate the standard deviation is below, where "d" represents the difference between the pair of duplicates, and "n" is the number of samples.

$$\text{Standard deviation } S = \sqrt{\frac{d^2}{2n}}$$

DISCUSSION

Patients with diabetes will continue to look for ways to cope with the high cost of managing the disease. One in every five dollars spent on healthcare in America is for the care of people diagnosed with diabetes. The average medical spending for a person with diabetes is about \$13,700, which is 2.3 times that of people without diabetes.

Due to cost, U.S. patients failed to fill 6.8% of the name-brand prescriptions their doctors requested in the 2008 fourth quarter, a 22% increase from the first quarter of 2007.

Actions like these come at the expense of patient health and add to the annual economic cost of the disease – \$245 billion spent on persons diagnosed with diabetes in 2012. This figure includes more than \$29 billion for diabetes supplies and agents, and \$44 billion for prescription medications. Indirect costs resulting from increased absenteeism, reduced productivity, disease-related unemployment disability, and loss of productive capacity due to early mortality add up to \$69 billion.

Blood glucose monitoring options don't have to be cost-prohibitive. Generic blood glucose test strips, also referred to as alternative test strips cleared by the US FDA for use with brand name meters-- such as UniStrip1™ Test Strips -- when used as intended, provide a high performance alternative that can save people with diabetes hundreds of dollars on test strips

alone (compared to name-brand strips).

The UniStrip1™ Test Strips have been subjected by the U.S. Food and Drug Administration to the same clinical testing, scrutiny and approval process as the name brands. The UniStrip1™ Test Strips have also been evaluated for accuracy and ease of use against the most recent and more restrictive ISO guidelines, and have met, or exceeded, the requirements for accuracy and precision.

Lower cost does not mean lower quality. This study demonstrates that the UniStrip1™ Test Strips deliver high levels of accuracy and precision when tested as intended. UNISTRIP TECHNOLOGIES products are offered for sale at price points significantly below those of name brands. People with diabetes can now use and rely on the UniStrip1™ Test Strips like they have become accustomed to doing with generic pharmaceuticals. The UniStrip1™ Test Strips are the low cost, high quality, alternative for blood glucose testing with the LifeScan OneTouch® Ultra® test strips with the OneTouch Ultra®, Ultra®2, UltraSmart®, and UltraMini® meters.

Study conducted by Neng-Chun Diabetes Clinical, Taiwan, ROC:
Yi-Ying Chu, I-Ching Wang and Dr. Neng-Chun Yu

Introduction and Discussion Sections by UniStrip Technologies,
LLC and Green Rhino Marketing, LLC

References:

http://www.cdc.gov/diabetes/pubs/pdf/ndfs_2011.pdf

IMS Health, National Sales Perspectives, Nov. 2011, National Prescription Audit, Dec. 2011

Wall Street Journal, April 8, 2009, Rockoff: "Many Drug Prescriptions are Now Going Unfilled"; with citations
to Wolter Kluwers Health

Economic Costs of Diabetes in the U.S. in 2012, Diabetes Care, published online ahead of print March 6,
2013; doi:10.2337/dc12-2625 Diabetes Statistics. Retrieved from <http://www.diabetes.org/diabetes-basics/diabetes-statistics> Chicago: "Diabetes Statistics." American Diabetes Association. Jan 15, 2014.

FOR ADDITIONAL INFORMATION OR QUESTIONS, CONTACT:

Unistrip Technologies, LLC
7224 Statesville Road, Suite A
Charlotte, NC 28269

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Technologies, LLC y Green Rhino Marketing, LLC